











BEGO compatibility overview 3D printing system components

Select your printer manufacturer

			
			
		<p>DLP printers in general</p>	







BEGO compatibility overview 3D printing system components



BEGO Varseo XS

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/¹ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/¹ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>											
 <p>VarseoWax Model Dental models</p>	<p>Ultrasonic bath: 5 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 5 min (Isopropanol 99%)</p>	<p>Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99% or Ethanol 96%)</p>									
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



BEGO Varseo/Varseo L/Varseo S

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/¹ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/¹ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>					<p>BEGO Otoflash: 2 x 500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>					
 <p>VarseoWax Model Dental models</p>												
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>	<p>Ultrasonic bath: 3 min + 2 min (Ethanol 96%)</p>					<p>BEGO Otoflash: 2 x 1.000 + 2 x 2.000 flashes</p>	<p>HiLite Power (Kulzer¹): 1 x 90 + 2 x 180 sec</p>					
 <p>VarseoWax Tray Individual impression trays</p>	<p>Ultrasonic bath: 3 min + 2 min (Ethanol 96%)</p>					<p>BEGO Otoflash: 2 x 2.000 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 180 sec</p>					

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.







BEGO compatibility overview 3D printing system components



Ackuretta¹ DENTIQ/Freeshape 120⁴

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>												
 <p>VarseoWax Model Dental models</p>	<p>Ultrasonic bath: 5 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 5 min (Isopropanol 99%)</p>	<p>Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99% or Ethanol 96%)</p>			<p>BEGO Otoflash: 2 × 2.000 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 180 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>Anycubic Wash & Cure Plus³: 2 × 20 min</p>			
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

⁴ Ackuretta Dentiq & Freeshape 120 can only be used with anodized aluminium build platform (for VarseoSmile Crown^{plus} and VarseoSmile Temp).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.

Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.

Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Ackuretta¹ SOL⁵

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox² Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>												
 <p>VarseoWax Model Dental models</p>												
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
² Both devices are identical in construction.
⁵ Ackuretta SOL can only be used with anodized aluminium build platform (Small, Medium and Large).
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Anycubic¹ Photon Mono X⁶

Cleaning

Light-curing

VarseoSmile Crown^{plus}
Permanent single crowns, inlays, onlays, and veneers

VarseoSmile Temp
Temporary crown and bridge restorations, inlays, onlays, and veneers

VarseoWax CAD/CAST
Burnout objects

VarseoWax Model
Dental models

VarseoWax Surgical Guide
Surgical guides and placement aids for implant prosthetics

VarseoWax Tray
Individual impression trays

<p>Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p>Formlabs¹ Form Wash: 5 min (Isopropanol 99 %)</p>	<p>Anycubic Wash & Cure Plus³: 8 min (Isopropanol 99 % or Ethanol 96 %)</p>

<p>Flash or LED light-curing device, e.g.:</p> <table border="0"> <tr> <td>BEGO Otoflash: 2 x 2.000 flashes</td> <td>HiLite Power (Kulzer¹): 2 x 180 sec</td> <td>Formlabs Form Cure: 2 x 20 min @ 60 °C</td> <td>Anycubic Wash & Cure Plus³: 2 x 20 min</td> </tr> </table>	BEGO Otoflash: 2 x 2.000 flashes	HiLite Power (Kulzer¹): 2 x 180 sec	Formlabs Form Cure: 2 x 20 min @ 60 °C	Anycubic Wash & Cure Plus³: 2 x 20 min
BEGO Otoflash: 2 x 2.000 flashes	HiLite Power (Kulzer¹): 2 x 180 sec	Formlabs Form Cure: 2 x 20 min @ 60 °C	Anycubic Wash & Cure Plus³: 2 x 20 min	

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
⁶ Compatibility applies to the design status to the serial no. P02123C0306508
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.









BEGO compatibility overview 3D printing system components



Asiga¹ Max UV/Max 405

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/¹ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: <i>See last page</i></p>	<p>SprintRay¹ ProWash/Dry: 4 min Reservoir 1 + 3 min Reservoir 2 + 3 min Trocknung + Absprühen (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/¹ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: <i>See last page</i></p>	<p>SprintRay ProWash/Dry: 4 min Reservoir 1 + 3 min Reservoir 2 + 3 min Trocknung + Absprühen (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>											
 <p>VarseoWax Model Dental models</p>	<p>Ultrasonic bath: 5 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 5 min (Isopropanol 99%)</p>	<p>Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99% or Ethanol 96%)</p>									
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Asiga¹ Pro 4K80⁷

Cleaning

Light-curing

VarseoSmile Crown^{plus}
Permanent single crowns, inlays, onlays, and veneers

VarseoSmile Temp
Temporary crown and bridge restorations, inlays, onlays, and veneers

VarseoWax CAD/CAST
Burnout objects

VarseoWax Model
Dental models

VarseoWax Surgical Guide
Surgical guides and placement aids for implant prosthetics

VarseoWax Tray
Individual impression trays

<p>Ultrasonic bath: 5 min (Isopropanol 99 % or Ethanol 96 %)</p>	<p>Formlabs¹ Form Wash: 5 min (Isopropanol 99%)</p>	<p>Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99 % or Ethanol 96 %)</p>

<p>Flash or LED light-curing device, e.g.:</p> <table border="0"> <tr> <td>BEGO Otoflash: 2 × 2.000 flashes</td> <td>HiLite Power (Kulzer¹): 2 × 180 sec</td> <td>Formlabs Form Cure: 2 × 20 min @ 60 °C</td> <td>Anycubic Wash & Cure Plus³: 2 × 20 min</td> </tr> </table>	BEGO Otoflash: 2 × 2.000 flashes	HiLite Power (Kulzer¹): 2 × 180 sec	Formlabs Form Cure: 2 × 20 min @ 60 °C	Anycubic Wash & Cure Plus³: 2 × 20 min
BEGO Otoflash: 2 × 2.000 flashes	HiLite Power (Kulzer¹): 2 × 180 sec	Formlabs Form Cure: 2 × 20 min @ 60 °C	Anycubic Wash & Cure Plus³: 2 × 20 min	

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
³ Compatibility applies to the design status to the serial no. W31126A0405446.
⁷ Printer must be operated in 4K mode.
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.



BEGO compatibility overview 3D printing system components

DEKEMA¹ trix print²

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>												
 <p>VarseoWax Model Dental models</p>												
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



EnvisionTEC¹ D4K Pro⁸

Cleaning

Light-curing

<p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>		Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)	Formlabs¹ Form Wash: 3 min (Isopropanol 99%)	Whip Mix¹ Veriwash + Veriwhirl/¹ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)	Manual cleaning with tensides: See last page	SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)	BEGO Otoflash: 2 x 1.500 flashes	HiLite Power (Kulzer¹): 2 x 90 sec	Formlabs Form Cure: 2 x 20 min @ 60 °C	SprintRay ProCure: 2 x 20 min @ 20 °C	SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown ^{plus} (60 sec – pause – 50 sec), Zone A	CUREbox¹ Plus: 2 x 20 min @ 30 °C	Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n
		Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)	Formlabs Form Wash: 3 min (Isopropanol 99%)	Whip Mix Veriwash + Veriwhirl/¹ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)	Manual cleaning with tensides: See last page	SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)	BEGO Otoflash: 2 x 1.500 flashes	HiLite Power (Kulzer¹): 2 x 90 sec	Formlabs Form Cure: 2 x 20 min @ 60 °C	SprintRay ProCure: 2 x 20 min @ 20 °C	SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A	CUREbox Plus: 2 x 20 min @ 30 °C	Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n
<p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>													
<p>VarseoWax CAD/CAST Burnout objects</p>													
<p>VarseoWax Model Dental models</p>													
<p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>													
<p>VarseoWax Tray Individual impression trays</p>													


¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.


² Both devices are identical in construction.







⁸ For the use of VarseoSmile Crown plus und VarseoSmile Temp please contact a distributor of EnvisionTEC GmbH.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components

Formlabs¹ Form 2⁹ (stainless steel build platform necessary) 

Formlabs Form 3B/3B+ (stainless steel build platform necessary) 

	Cleaning	Light-curing
 VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers		
 VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers	Formlabs Form Wash: 3 min (Isopropanol 99%)	Formlabs Form Cure: 2 x 20 min @ 60 °C
 VarseoWax CAD/CAST Burnout objects		
 VarseoWax Model Dental models		
 VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics		
 VarseoWax Tray Individual impression trays		

	Cleaning	Light-curing
	Formlabs Form Wash: 3 min (Isopropanol 99%)	Formlabs Form Cure: 2 x 20 min @ 60 °C
	Formlabs Form Wash: 3 min (Isopropanol 99%)	Formlabs Form Cure: 2 x 20 min @ 60 °C

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.
⁹ VarseoSmile Crown^{plus} is distributed by Formlabs as Permanent Crown + VarseoSmile Temp as Temporary CB.
 For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.
 Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.
 Not all products and services shown are available in all countries.







BEGO compatibility overview 3D printing system components



Microlay¹ Versus 385

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: <i>See last page</i></p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: <i>See last page</i></p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>												
 <p>VarseoWax Model Dental models</p>	<p>Ultrasonic bath: 5 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 5 min (Isopropanol 99%)</p>	<p>Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99% or Ethanol 96%)</p>			<p>BEGO Otoflash: 2 × 2.000 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 180 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>Anycubic Wash & Cure Plus³: 2 × 20 min</p>			
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



SprintRay¹ Pro 95/Pro 95 S/Pro 55/Pro 55 S

Cleaning

Light-curing

<p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix¹ Veriwash + Veriwhirl/ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
<p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 × 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 × 90 sec</p>	<p>Formlabs Form Cure: 2 × 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 × 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 × 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 × 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
<p>VarseoWax CAD/CAST Burnout objects</p>												
<p>VarseoWax Model Dental models</p>												
<p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
<p>VarseoWax Tray Individual impression trays</p>												

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

BEGO compatibility overview 3D printing system components



Whip Mix¹ VeriBuild/VeriEko¹⁰

Cleaning

Light-curing

 <p>VarseoSmile Crown^{plus} Permanent single crowns, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs¹ Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/¹ Ackuretta¹ Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay¹ ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Crown^{plus} (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox¹ Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoSmile Temp Temporary crown and bridge restorations, inlays, onlays, and veneers</p>	<p>Ultrasonic bath: 3 min + 2 min (Isopropanol 99% or Ethanol 96%)</p>	<p>Formlabs Form Wash: 3 min (Isopropanol 99%)</p>	<p>Whip Mix Veriwash + Veriwhirl/¹ Ackuretta Cleani²: 3 min + 3 min (Isopropanol 99%)</p>	<p>Manual cleaning with tensides: See last page</p>	<p>SprintRay ProWash/Dry: 4 min reservoir 1 + 3 min reservoir 2 + 3 min drying + spraying off (Isopropanol 99%)</p>	<p>BEGO Otoflash: 2 x 1.500 flashes</p>	<p>HiLite Power (Kulzer¹): 2 x 90 sec</p>	<p>Formlabs Form Cure: 2 x 20 min @ 60 °C</p>	<p>SprintRay ProCure: 2 x 20 min @ 20 °C</p>	<p>SprintRay ProCure 2: Preprogrammed resin profile VarseoSmile Temp (60 sec – pause – 50 sec), Zone A</p>	<p>CUREbox Plus: 2 x 20 min @ 30 °C</p>	<p>Ackuretta CURIE: 2 x 2,5 min Exposure parameters: P13 D8 T2.30B0n</p>
 <p>VarseoWax CAD/CAST Burnout objects</p>												
 <p>VarseoWax Model Dental models</p>												
 <p>VarseoWax Surgical Guide Surgical guides and placement aids for implant prosthetics</p>												
 <p>VarseoWax Tray Individual impression trays</p>												

NEW!
Also available in Bleach

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

² Both devices are identical in construction.

¹⁰ Whip Mix VeriBuild & VeriEKO can only be used with anodized aluminium build platform (VeriEKO anodized aluminium build platform in Small, Medium and Large).

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins.

Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality.

Not all products and services shown are available in all countries.



BEGO compatibility overview 3D printing system components

DLP printers in general (with 385 – 405 nm wave length)

Cleaning

Light-curing



VarseoSmile Crown plus
Permanent single crowns, inlays, onlays, and veneers



VarseoSmile Temp
Temporary crown and bridge restorations, inlays, onlays, and veneers



VarseoWax CAD/CAST
Burnout objects



VarseoWax Model
Dental models



VarseoWax Surgical Guide
Surgical guides and placement aids for implant prosthetics



VarseoWax Tray
Individual impression trays

--	--	--

--	--	--

--	--	--

Ultrasonic bath: 5 min (Isopropanol 99% or Ethanol 96%)	Formlabs¹ Form Wash: 5 min (Isopropanol 99%)	Anycubic¹ Wash & Cure Plus³: 8 min (Isopropanol 99% or Ethanol 96%)
---	--	---

--	--	--

--	--	--

--

--

--

Flash or LED light-curing device, e.g.:
BEGO Otoflash: 2 × 2.000 flashes
HiLite Power (Kulzer¹): 2 × 180 sec
Formlabs Form Cure: 2 × 20 min @ 60 °C
Anycubic Wash & Cure Plus³: 2 × 20 min

--

--

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group.

³ Compatibility applies to the design status to the serial no. W31126A0405446.

For detailed information on the production workflow, please refer to the instructions for use of the respective VarseoSmile or VarseoWax resins. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.

Manual cleaning with Tensides

This cleaning method is valid for VarseoSmile Crown^{plus} and VarseoSmile Temp


Necessary tools, equipment and materials

- InovaPrint wash (hp-dent¹) general purpose 3D print cleaner
- Tap water
- 1-propanol (70 Vol.-%)
- Toothbrush
- Brush
- Instrument for holding the printed objects (e.g. artery clamp)
- Absorbent pad (e.g. paper towels) or tub
- 2 cups for 1-propanol (70 Vol.-%)
- 1 cup for cleaning solution
- Compressed air with trigger/splash guard
- Personal protective equipment: protective gloves and goggles



Cleaning process

Carefully remove the printed object from the platform with the help of a spatula.



Clean the printed object by using a toothbrush and/or brush in two steps:

Step 1:

Pre-wash with reusable cleaning solution: 5 % InovaPrint wash (hp-dent) + 95 % tap water.

- 1.1 Vigorously swirl object in cleaning solution for 15s using pliers or an arterial clamp. Ensure not to damage the printed object.



- 1.2 Remove excess resin by using a toothbrush and cleaning solution. In addition, a brush can be used to clean the inside of crowns. Shortly swirl object in cleaning solution as needed.



Note: Clean toothbrush and brush regularly with 1-Propanol when covered with resin. Dry before reusing toothbrush or brush again for cleaning the printed object.

- 1.3 Use compressed air under a fume hood with splash guard to remove the cleaning solution from the printed object's surface.



- 1.4 Repeat step 1.2 and 1.3 until just a thin layer of resin remains on the surface of the print.

Step 2:

Finish with fresh 1-Propanol (70 Vol.-%).

- 2.1 Vigorously swirl printed object in 1-Propanol for 5s and dry object immediately with compressed air.



- 2.2 Check for remaining resin (shiny spots). If surface of printed object is matte, remove support structure as described in the following step (step 2.3). Shiny spots can be removed by quickly brushing the surface with a 1-Propanol soaked brush. Immediately dry object with compressed air.

- 2.3 Remove the support structure with the help of a cutting wheel or side cutters. Ensure not to deform the printed object.



- 2.4 Quickly clean occlusal surface with a 1-Propanol soaked brush to remove any excess resin. Immediately dry printed object with compressed air.



Note: Contact time of printed object with 1-Propanol needs to be reduced to a minimum to avoid the formation of white spots on the surface.

- 2.5 Post-cure the dental objects by using a validated post-curing unit.

¹ This symbol is a commercial designation/registered trademark of a company which is not part of the BEGO company group. Pictures and illustrations are exemplary. Colors, symbols, design and information on the labels and/or packaging shown may differ from reality. Not all products and services shown are available in all countries.